

AGRI-KNOWS

Knowledge Transfer in Agriculture as Added Value in Environment Protection

Javni razpis/ AGRI - KNOWS/ Prenos znanja v kmetijstvu kot dodana vrednost pri zaščiti okolja

Javni razpis/ AGRI - KNOWS/Trasferimento delle conoscenze in agricoltura come valore aggiunto per la tutela dell'ambiente

Public procurement/ AGRI - KNOWS/ Knowledge transfer in agriculture as added value to environment protection

Karmen Goljevšček Čargo

Translated by Klara Levstek

Šempeter pri Gorici, 28.2.2014



2007-2013

cooperazione territoriale europea
programma per la cooperazione
transfrontaliera

Italia-Slovenia

evropsko teritorialno sodelovanje
program čezmejnega sodelovanja

Slovenija-Italija



Investiamo nel
vostro futuro!

Naložba v vašo
prihodnost!

www.ita-slo.eu

Progetto cofinanziato dal Fondo europeo di
sviluppo regionale


Projekt sofinancira Evropski sklad
za regionalni razvoj

THE PURPOSE OF THE RESEARCH

The purpose of the research was to see how susceptible the microorganisms are to penicillin. For this purpose we used antibiogram.



PROCEDURE




ANTIBIOGRAM



PREPARING SOIL AND WATER
SAMPLES



INOCULATION OF WATER AND SOIL
MICROORGANISMS ON SOLID MEDIA



COUNTING COLONIES ON
INDIVIDUAL GROWTH MEDIA

I. ANTIBIOGRAM

- *Antibiogram is a method which is used for defining bacterial susceptibility to antibiotics or their resistance to them.*
- *The substances used in that method are chemotherapeutics, which destroy microorganisms. Antibiotics are one type of them.*



We distinguish 3 antibiograms:

ANTIBIOGRAM BASED ON DISC DIFFUSION METHOD

(Antibiotic discs or strips are applied on a solid culture medium.)

ANTIBIOGRAM BASED ON DILUTION METHOD (Bacterial growth is observed in a series of liquid media with different antibiotic concentrations.)

CUMULATIVE ANTIBIOGRAM (Different bacteria are inoculated on a bacterial growth medium, which contains an antibiotic.)




Procedure:

-In a sterile Erlenmeyer flask we prepared a soil suspension (10 g soil + 90 ml sterile physiological solution).

-We shook the suspension on a chemical shaker for 1 hour.

-In test tubes, we prepared 6 different penicillin solutions: from 10^{-1} to 10^{-6} .









1 ml penicillin
+9 ml sterile physiological solution
 10^{-1}

- *With a pipette, we applied 0.1 ml soil suspension on each of the 6 media.*
- *Using a spatula, we spread the soil sample evenly over the whole surface of a petri dish. Afterwards we dried it by the fire for 5 minutes.*
- *Using tweezers, we put into each petri dish a disc, which had been soaked in one of the different penicillin solutions before.*





Results:



CHOOSING A GROWTH MEDIUM

GROWTH MEDIA are substrates for growing microorganisms in lab conditions.

We distinguish different culture media:


- *COMPLEX MEDIA: They are rich complex media where most heterotrophic bacteria can grow.*
- *SELECTIVE MEDIA: They contain ingredients that encourage the growth of specific bacteria and inhibit the growth of others.*




- *In our research we used two types of growth media:*

a) nutrient broth (PKE)

*b) MacConkey Agar
(a selective growth medium for coliform bacteria)*





PRENOS ZNANJA V KMETIJSTVU KOT DODANA VREDNOST V VAROVANJU OKOLJA

Javni razpis št. 03/2009/ Bando publico n. 03/2009: AGRI - KNOWS

Karmen G.Čargo/ Šempeter pri Gorici, 28.2.2014




Materials:

- a soil or water sample
- A sterile medium
- An inoculation loop
- A lab burner



How are microorganisms inoculated on a medium?

All microbiological tests must be done aseptically.



This means that we deal with all contaminants and bacterial cultures in such a way that there is no invasion of undesired microorganisms. At the same time we take care not to spread any microorganisms.

Glassware, objects and utensils as well as growth media have to be **sterile**.



We use the flame from a gas burner.

We flame sterilize the inoculation loop and glassware.



II. INOCULATION OF MIKROORGANISMS ON SOLID MEDIA (WATER)

This experiment was done out to see how much penicillin is leached into water.



1.

- Making columns from plastic tubes, filling them with soil

2.

- Polluting the columns with penicillin and pouring rainwater over

3.

- Inoculating microorganisms from rainwater on a solid medium, counting colonies, analyzing all samples

4.

- Inoculating microorganisms from soil on a solid medium, counting colonies, analyzing all samples



1. Making the columns

PROCEDURE


- We took six columns (plastic tubes): 10 cm in diameter, 3 mm thick and 50 cm high. All the columns had a cone shape on the bottom.
- At the bottom of the column, we place a piece of plastic net and pour Over 40 ml coarse sand. Then we add another piece of plastic net and cover it with 50 ml fine sand.



- Afterwards we gradually filled in the columns with 4 kg soil, which we had taken from our school green house. After each kilo we stopped and compressed the soil. Then we continued.
- Then we attached a plastic hose to the conical bottom of the columns. The hose was about 150 cm long.
- Finally, we poured water over the columns. We had to make sure that the water height in the column and in the hose was 60 cm. In that way we reached the same pressure. After two days we released the water out of the hose.



2. Polluting the columns with penicillin



3. Preparing water samples

- Into each column, we poured 400 ml rainwater. This rainwater went through the soil and slowly flowed out through the hose.



- We collected that water in a measuring cup. Then we analyzed it.
- We repeated the procedure 4 times.





PRENOS ZNANJA V KMETIJSTVU KOT DODANA VREDNOST V VAROVANJU OKOLJA

Javni razpis št. 03/2009/ Bando publico n. 03/2009: AGRI - KNOWS
Karmen G.Čargo/ Šempeter pri Gorici, 28.2.2014



4. Inoculating microorganisms on the medium

- We used the flame from a gas burner. We flame sterilized the inoculation loop and glassware.
- We regularly cleaned the worktop with disinfectants.
- We inoculated microorganisms from rainwater on both growth media, using a flame sanitized inoculation loop.




- We took one drop of water aseptically on the growth medium.

- We spread the water drop over the petri dish with wavy movements.

- Afterwards we flame sterilized the inoculation loop until it started glowing.

- Additionally, we applied another water drop on the growth medium and made a wavy line, which was perpendicular to the first one.

We got a twisted pattern.



5. The growth of microorganisms:

- We incubated the growth medium at 28 °C.



6. Results:



III. INOCULATION OF MICROORGANISMS ON SOLID GROWTH MEDIA (SOIL)

This experiment was carried out in order to see how many antibiotics stayed in soil.



1. Preparing a soil sample

- In an Erlenmeyer flask we weighed 10 g soil and poured over 90 ml sterile physiological solution.
- We put the Erlenmeyer flask on a chemical shaker at speed 2 for one hour.




2. Inoculating microorganisms on the medium

-From that suspension we inoculated microorganisms on both growth media, using a flame sanitized inoculation loop.



- We spread the drop over the petri dish with wavy movements.
- Afterwards we flame sterilized the inoculation loop until it started glowing.
- Additionally, we applied another water drop on the growth medium and make a wavy line, which was perpendicular to the first one.




3. The growth of microorganisms:

- We incubated the growth medium at 28 °C.



4. Counting colonies

- The number of colonies which grew on each medium



IV. RESULTS

NUM BER	DATE OF WATER ANALYSIS	DATE OF SOIL ANALYSIS
I.	18/4/2013	18/4/2013
II.	26/4/2013	
III.	23/5/2013	
IV.	28/5/2013	28/5/2013



	1ST COLUMN		2ND COULMN		3RD COULMN		4TH COULMN		5TH COULMN		6TH COULMN	
	PKE	KOLI	PKE	KOLI	PKE	KOLI	PKE	KOLI	PKE	KOLI	PKE	KOLI
I.	243	46	195	80	101	14	108	10	167	9	131	9
II.	200	unco unta ble	153	141	135	26	141	19	157	8	unco unta ble	8
III.	71	87	81	102	155	60	176	45	160	12	131	12
IV.	uncou ntable	174	unco untab le	uncou ntable	unco untab le	176	uncou ntable	unco unta ble	unco unta ble	22	unco unta ble	22



AGRI-KNOWS

Prenos znanja v kmetijstvu
kot dodana vrednost
v varovanju okolja

Šempeter pri Gorici, 28.2.2014

Hvala za vašo pozornost!

Grazie per l'attenzione!

Thank you for your attention.

karmen.goljevscek@guest.arnes.si



Projekt [AKRONIM/NASLOV PROJEKTA] je sofinanciran v okviru Programa čezmejnega sodelovanja Slovenija-Italija 2007-2013 iz sredstev Evropskega sklada za regionalni razvoj in nacionalnih sredstev

Progetto [inserire l'acronime/titolo del progetto] finanziato nell'ambito del Programma per la Cooperazione Transfrontaliera Italia-Slovenia 2007-2013, dal Fondo europeo di sviluppo regionale e dai fondi nazionali.



REPUBLIKA SLOVENIJA
MINISTRSTVO ZA GOSPODARSKI
RAZVOJ IN TEHNOLOGIJO



Ministero dell'Economia
e delle Finanze



2007-2013

cooperazione territoriale europea
programma per la cooperazione
transfrontaliera

Italia-Slovenia

evropsko teritorialno sodelovanje
program čezmejnega sodelovanja

Slovenija-Italija



Investiamo nel
vostro futuro!

Naložba v vašo
prihodnost!

www.ita-slo.eu

Progetto cofinanziato dal Fondo europeo di
sviluppo regionale

Projekt sofinancira Evropski sklad
za regionalni razvoj